

# 新型コロナウイルスPCR検査 | 検体採取マニュアル


Setolabo 衛生検査所



検体受付フォーム

QRコードを読み取ってください。

## 1. 検査キットを受け取り、中身を確認する



検査キット

XX0000


検査番号シール

田中太郎

090-0000-0000

個人情報シール

2枚貼ってあるか  
確認する。



検査番号カード  
(お控えの紙)

検査番号シールと  
カードの検査番号が  
合っているか確認する。

カードは持って帰り  
**検体受付フォーム**に  
必ず登録をする。  
**(QRコード)**


## 2. 唾液をチューブの中に入れてみましょう



- ① 唾液採取用チューブの蓋を開けます。
- ② 喉にウイルスが溜まりやすいので、喉元を軽くマッサージしてから唾液を採取してください。

唾液を容器の底から  
1cm入れてください。  
ふたを最後まで  
しっかりと締めます。

※ 「ふた」や「チャック付きバック」  
は、テープでぐるぐる巻きにしない  
てください！  
中身が漏れないように閉まっ  
ているか、必ず確認してください。



唾液採取前に歯磨きをして、**30分**時間を空けてから採取してください。

## 3. 検体採取後、窓口スタッフに提出する


OK

正しいセットの見本



NG

悪い例




※ 下記のようにしないでください！

- ふたがしっかり閉まっていない。
- ふたをテープでぐるぐる巻きにしている。
- 唾液が1cm入っていない。
- チャック付きバックをテープでぐるぐる巻きにとめている。




## 1. Receive the test kit and check it. / 收到测试包并检查内容。



Test kit / 检查套件

**XX0000**


Test number label / 检验编号标签

**田中太郎**  
**090-0000-0000**Personal information label  
/ 个人信息标签Check if two label are pasted.  
检查是否粘贴了两张标签。Test number card  
检验编号卡Check that the Test number label  
matches the Test number card.  
确保 [ 检验编号标签 ] 和  
[ 检验编号卡 ] 的检验编号匹配。Take the card home and  
be sure to register it on the  
**Specimen reception form.**  
**(QR code)**  
确保将卡带回家并在  
**标本接收表**上进行注册。  
**(二维码)**



## 2. Put saliva in the tube. / 让我们把口水放进管子里。



- ① Open the lid of the tube. / 打开管的盖子。
- ② Viruses tend to collect in the throat, so massage the throat lightly before collecting saliva. 病毒往往会聚集在喉咙中，因此在收集唾液之前请轻轻按摩一下喉咙。

Put saliva 1 cm from the bottom of  
the container. Tighten the lid tightly.  
在容器底部放 1 厘米的唾液。拧紧盖子。**Do not wrap the "lid" or "zipper bag" with tape!** Be sure to check that it is closed so that the contents do not leak.**不要用胶带包裹“盖子”或“拉链袋”！** 确保检查它是否已关闭，以防止内容物泄漏。After brushing your teeth, Please be saliva taken from spaced **30 minutes** or more time.  
刷牙和收集唾液等待至少 **30 分钟**。

## 3. After specimen collection, submit it to staff. / 收集样品后，将其提交给工作人员。

**OK****Correct set**  
一组正确的样品**NG****Bad example / 不好的例子****※ Please do not below! / 请不要做以下!**

- The lid is not closed tightly. / 盖子未关紧。
- The lid is taped. / 盖上胶带。
- Saliva is less than 1 cm. / 唾液少于 1 厘米。
- The zipper bag is wrapped around with tape. / 拉链袋用胶带包裹。